

NQ09

Portable Heated Incubator



Part #	Description	Function	Temp Range	Interior (Metric)	Exterior (Metric)	Weight (Metric)	NQRB Duration
NQ09	Portable Incubator	Heating	Ambient+5C to 40C	13 cm (H) 27 cm (W) 19 cm (D)	27 cm (H) 33 cm (W) 24 cm (D)	3.6 kg	8 hours

*NQRB duration only applies to incubators equipped with rechargeable battery packs.

Durations are measured after the chamber has been preconditioned.